REDAX®

IT

VALVOLA DI HEIMLICH

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DESCRIZIONE

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EN

DESCRIPTION

The Heimlich non-return valve is a one-way drainage device which prevents the air reflux towards the patient's chest.

Two versions of the device are available:

Standard Heimlich valve (fig. 1)

This device is composed of the following parts:

1. A flexible body (1) fitted with universal tapered connectors (2) at both ends, so as to be connected to drainage tubes and/or collection devices.
2. A lip valve in flexible material (3), connected to one of the two connectors.
3. An additional connector (4) for the connection to the drainage tubes provided with female luer-lock connector.

Heimlich valve with bag (fig. 2)

This device is composed of the following parts:

1. A standard Heimlich valve (as described in the previous paragraph), the end of which is connected to a collection bag (5), by means rubber connection tube (6).
2. The collection bag has 2000 ml capacity and it is provided with an inlet one-way valve (7).
3. Furthermore, the bag is provided with a filter (8) that allows to evacuate air from the patient outwards.

RECOMMENDATIONS FOR USE

The Heimlich valve is primarily indicated in drainage of the chest cavity by gravity in order to allow the evacuation of liquids and/or air, as well as preventing the backflow towards the patient. This way, the normal gradients of intrathoracic pressure can be restored and it can facilitate a complete pulmonary expansion, by restoring the normal respiratory dynamics. Both versions, with and without bag, are especially indicated in emergency conditions.

Note: This device can only be used in safe conditions in the field of applications and using the procedures indicated on this operating instruction sheet, for the specific type of product. The manufacturer declines all liability relating to improper use or use different to that indicated.

Fig. 1


Fig. 2



Fig. 3

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