

MANUFACTURER DECLARATION

To whom it may concern,

We OSATU, S.coop. (trademark Bexen Cardio) manufacturer of medical devices, with address in Edificio Zearekobuelta, Subida de Areitio nº 5, 48260 Ermua (Vizcaya) SPAIN,

Hereby declare that:

The REANIBEX 100 performs automatic self-tests when switched on, during operation and periodically when the equipment is powered or switched off. Regular routine maintenance is recommended to complement the self-tests performed by the equipment, which minimally includes:

- Make sure that the REANIBEX 100 is clean and that neither the equipment nor its accessories are visibly damaged.
- After each use, ensure that all the necessary supplies, accessories and spare parts are available for the next use.
- With the cover closed, make sure that the status indicator LED flashes every 10 seconds. If not, open the cover and check if the pads are correct. Perform a Manual test to check the equipment status.

Signed in Ermua, October 17, 2024


Bexen CARDIO
Osatu S.Coop. - F48113161
Edificio Zearekobuelta
Subida de Areitio 5, 48260 Ermua (Spain)
www.bexencardio.com
Oxel Etxebarria
Managing Director